



INSTRUCTION GUIDE

Microbiome Explorer[™]

Comprehensive



We're happy you're here!
At Microba, we expect better –
better for you and better for your
health. We're with you as you take
the positive steps towards clearer
insights into your microbiome.

For questions or further
information on what
happens next, visit our
support page or call us
on **+44 808 501 2748**.



What to know two weeks ahead of taking your sample

When to collect and return your sample

Only collect your sample on a Monday, Tuesday, or Wednesday. Return it to your local post office before 1pm on the same day.*

Continue your usual routine

Keep your regular exercise, diet, and medications unless instructed differently by your healthcare provider.

What kind of sample to provide

Take the sample from a normal bowel movement that reflects how things usually are for you.

If you menstruate

Avoid collecting your sample during your period or within three days afterward.

If you're on medication

Please note that antibiotics, antimicrobials and other medications can impact test results. Please refer to our FAQs at the back of this booklet for more information, or consult your healthcare provider.

*Post office express post packages cut-off time



48 hour Preparation

When receiving kit



- Ensure the kit is addressed to you. Your sample is identified by the QR code on the swab base, the top of the pot, and your completed collection stickers.
- Unpack your kit and check you have all the components.
- Read and understand all instructions.
- Freeze your ice pack for sample return in the box.

Two days before taking your sample



- Avoid drinking alcohol.
- Avoid taking NSAIDs (e.g. ibuprofen, aspirin) unless directed by a healthcare professional.



24 hour Preparation

One day before taking your sample

- Freeze your ice pack for sample return in the box.
- Be ready to take your sample from your first bowel movement of the day.
- Allow 30 minutes to complete the sampling process.


Day of sample

- Collect your sample on a **Sunday, Monday, or Tuesday**. Drop it off before 1 pm the same day, or on Monday morning if you collect it on Sunday.
- For women, avoid collecting during or three days after menstruation.
- Ensure you have read and understood all instructions. If you need help, call Customer Services at **+44 808 501 2748**.




Make sure you have all of the below items included in your Microbiome Explorer™ kit

The below items are included in your sampling kit:



1 Instruction booklet



2 Checklist



3 Collection paper



4 Device 1



5 Device 2



6 Device 3



7 Ice pack



Before starting your test, ensure the kit is addressed to you and contains all listed items. If anything is missing or damaged, contact Microba Customer Service at hello.uk@microba.com or call +44 808 501 2748.




8 Collection stickers x 2



9 Pouch A




10 Pouch B



11 Kit box



12 Reply-paid envelope



13 Thermal pouch



STEP 1

Prepare kit items to take your sample

When you are ready to collect a sample, take the collection kit to the toilet. You will need:

- a. Instruction booklet
- b. Collection paper
- c. Capped vial (Device 1)
- d. Pot (Device 2)
- e. Swab (Device 3)





STEP 2

Prepare a collection paper on the toilet

Open and unfold the collection paper provided and slip it over the toilet seat.

Alternatively, you may line a clean disposable container with the collection tissue and place the container into the toilet bowl after urination. This collection method is recommended if you have loose bowel movements.



STEP 3

Provide faecal sample onto the collection paper

Ensure your faecal sample remains uncontaminated from urine, toilet water and other contaminants. You may need to urinate before opening the collection paper accessory and starting the faecal collection process.

In case you have a problem with your first collection attempt, use a disposable container lined with toilet paper.





STEP 4

Use Device 1 to collect your capped-vial sample

Write your full legal name, date of birth, and the date of collection on the label using a black or blue pen.

Twist the cap to open the sample bottle. Scrape the stick along the stool until the grooves are covered. You only need a small amount for the test. Please do not add extra. Place the stick back in the bottle and twist until secure. **Do not reopen the bottle after sealing it.**




STEP 5

Use Device 2 to collect your pot sample

Open Device 2 and use the scoop inside the pot lid to collect the faecal material. Collect material from three to five different sections of the sample to fill the pot with at least 1 tablespoon of sample. Make sure you don't fill the pot over halfway.


If blood or mucus is present in your stool, please avoid collecting from those areas.



STEP 6

Open the tube by twisting and pulling


Open Device 3 by twisting and pulling with a firm grip. It is important that the swab tip does not come into contact with your hand or any other surfaces, as this can compromise the sample.



Take a small amount of the stool from Device 2 using the swab. Only a small smear is needed with no lumps captured on the swab.

Immediately place the swab back in the tube without touching any other surfaces. Firmly close the tube until you hear it click.

Tear the edges of the collection tissue and drop into the toilet to flush.





Lightly touch the end of the swab to pick up a small smear – only a small discolouration is needed.

Avoid providing too much or too little sample on the swab, or our lab may not be able to process your sample.



Too much
sample



Too little
sample



STEP 7

Label your sample

Use the collection stickers provided to label Device 2 and Device 3 with your full legal name, date of birth, time (AM or PM), and date of collection.




STEP 8

Complete your Patient Consent Form

To ensure we can process your sample, please complete the patient consent form and return it along with your labelled sample.

We will not be able to process your samples without this information.




Consent For Sample Processing

I hereby grant Invivo Clinical Ltd and Microba Life Sciences Ltd the authority to process my biosample and personal data for the purposes of preventative and occupational medicine and research.

I am aware that I may withdraw my consent at any time by informing Invivo Clinical Ltd and Microba Life Sciences Ltd.

I also consent to the disclosure of the results to the referring healthcare provider, a member of their clinic, a clinician of my choice, or directly to myself. I can confirm that Invivo Clinical Ltd or Microba Life Sciences Ltd have not asked me to discontinue treatment of care from my primary care provider.


Print Name	John Doe
Date of Birth DD/MM/YY	5 / 3 / 1984
Signature	
Date Signed DD/MM/YY	15 / 2 / 2025



STEP 9

Place Device 1 and Device 3 in Pouch A


Place Device 1 and Device 3 in Pouch A provided and press the pouch closed to seal.



STEP 10

Place Device 2 in Pouch B


Place Device 2 in Pouch B provided. Do not discard the absorbent paper. Insert and press the pouch closed to seal.



STEP 11

Refrigerate Pouch A and Pouch B for a minimum of 30 minutes

Place Pouch A and Pouch B in the fridge for a minimum of 30 minutes or until ready to post that same day. Keep the ice pack in the freezer so it's ready to go when you send your sample.





STEP 12

Remove pouch B from the fridge and activate all temperatures sensors

Remove pouch B from the fridge when ready to post, ensuring the sample has been refrigerated for at least 30 minutes.



Activate the temperature sensors attached to the pouch in a room that is 30°C or below. To activate, fold up and pull the indicator activation tab until the tab and barrier film has been removed.



STEP 13

Pack the samples for return

Pack Pouch A, Pouch B, and the ice pack into the thermal pouch. Then place the thermal pouch inside the kit box (with the insert removed). Finally, place the kit box into the prepaid envelope.



STEP 14

Post over the counter at Post Office

Return to our laboratory by taking to your nearest Post Office and returning over the counter.



Frequently Asked Questions

Should I discontinue probiotics or supplements before sampling?

We recommend maintaining your normal supplement and probiotic regime for two weeks prior to sampling, unless you have been advised to discontinue them by your healthcare professional. Always consult with your healthcare professional before making any changes to prescribed probiotics or supplements and do not discontinue any medications without their advice. Please collect a sample from a bowel movement that is typical for you.

Will taking antibiotics and antimicrobial supplements prior to sampling change my results?

Antibiotics and antimicrobial supplements can disrupt the microorganisms in your gut microbiome. We recommend waiting four weeks or longer after ceasing your supplement or antibiotic course before taking your sample.

Can I provide a sample with visible blood present?

Yes, you can provide a sample with visible blood, though it's preferred to wait until bleeding resolves. Visible blood may result in a positive faecal occult blood reading. Confirm with your practitioner to ensure proper interpretation.



Should I stick to my regular diet and lifestyle habits before sampling?

If you would like to see what your 'normal' gut microbiome looks like, stick to your normal habits for at least two weeks prior to sampling. The microorganisms inhabiting your gut are quite stable over time (in the absence of any major disruptions) and it is only their abundance that will change based on your diet and other habits.

Can I provide a liquid stool sample?

Yes, you can still send a liquid stool sample for analysis. However, a liquid stool sample may impact result accuracy which will be noted on your report.

Can I provide a sample if I have been taking nonsteroidal anti-inflammatory drugs (NSAIDs)?

It is better to take a sample when you have not been taking NSAIDs in the few days prior. You can still send a sample if you have been taking NSAIDs, but this will impact the level of calprotectin reported. Please advise your healthcare professional so this can be considered in the report interpretation.

For answers on more Frequently Asked Questions visit our website www.microba.com or contact the Microba Customer Service team on +44 808 501 2748 for more information.



Symbol Glossary



In vitro diagnostic
medical device



Batch Code



Consult instructions for use



Do not re-use



Keep dry



Authorized representative
in the United Kingdom



Catalogue number



Manufacturer



Storage conditions



Keep away from sunlight



PRODUCT STORAGE

The collection kit must be stored in its original package at a temperature between 5°C and 30 °C until the time of use.

If incorrectly stored, its effectiveness will be compromised.

LIMITATIONS

Microbiome Explorer Sample Collection Kit is not intended for users under 18 years of age.

INTENDED USE

The Microba Microbiome Explorer Specimen Collection Kit is intended for lay adult users to self-collect a stool specimen at home. The kit is a single use in vitro diagnostic (IVD) designed to provide containment, packaging and method of transport for the specimen to a designated laboratory where analysis is performed.

WARNINGS AND PRECAUTIONS:

1. For in vitro diagnostic use.
2. Do not sample faeces from the toilet, faeces mixed with urine, menstrual blood or water.
3. Caution is recommended in using antacids, barium, bismuth, anti-diarrheal medication, antibiotics, histamine, nonsteroidal anti-inflammatory drug or laxatives prior to collection of the specimen.
4. Do not leave the device unattended.
5. Not suitable for rectal sampling.
6. Do not use after the expiry date printed on either device 1 or 3.
7. This product is for single use only. Reuse may cause inaccurate results.
8. Not suitable for any other application than intended use.





Microbiome Explorer is only available for purchase through a healthcare professional.

The faecal occult blood, polymerase chain reaction (PCR) and enzyme-linked immunosorbent assays (ELISA) used in Microbiome Explorer are diagnostic and are approved for clinical use. The faeces pH assay used in Microbiome Explorer is for research use, only and not to be used as a basis for diagnosis. The metagenomic assays used in Microbiome Explorer are to determine the microbiome populations and associated functional pathways in a faecal sample. The application is for research use only and is not to be used as a basis for diagnosis. Learn more about the journey we are on to validate this gold-standard technology for clinical diagnosis and application at microba.com.



Microba Pty Ltd
Level 10, 324 Queen Street, Brisbane, QLD, 4000, Australia



Invivo Clinical Ltd

The Coach House, 3 Lewiston Mill Toadsmoor Road, Stroud, GL5 2TE



MEK-003