

Deep Dive: Measuring the microbiome as an ecosystem

A healthy gut microbiome is a balanced one. Evaluating that balance cannot be reduced to a list of well-known species. It reflects the overall state of the ecosystem shaped by diversity, relative abundance, and the functional capacity of the microbial community.

Many testing approaches still measure only fragments of this system – selected pathogens, limited taxa, or small gene regions and attempt to infer balance from that partial view.

If the clinical question is ecosystem-level balance, the measurement must reflect the complexity of the ecosystem.



Partial testing can't reliably answer critical clinical questions

In clinical practice, microbiome testing is rarely about detection alone. The goal is interpretation: understanding how microbial ecology may be influencing:



Interpretation requires context. A partial view of the microbiome may detect something, but it cannot reliably answer the more clinically meaningful questions:

How does this organism's abundance compare to the overall ecosystem?

Is the ecosystem resilient or vulnerable to instability?

Is functional capacity reduced?

Is it contributing to an anti-inflammatory or pro-inflammatory state?

Why Targeted Testing Misses the Bigger Picture

Because microbiome balance emerges from the whole system, targeted tests are inherently limited in what they can explain. Healthy gut microbial communities are characterised by high taxa diversity, high microbial gene richness, and a stable functional core, properties that cannot be captured by targeted detection methods alone.¹



Refer a Patient

Let's deliver better, together.

[Order a test](#)

Testing technologies vary in what they can measure and reveal

Microbiome testing has progressed through several technological stages. Each has improved capability, but not all are designed to assess ecosystem-level balance. Understanding these differences helps clinicians choose appropriate tools and interpret results correctly.

Culture-Based Testing

Useful for pathogen detection, not for ecosystem assessment.

Culture remains valuable for identifying certain infectious organisms and guiding antimicrobial therapy²

However:

- Many gut organisms cannot be cultured using standard laboratory techniques.
- Culture favours organisms that grow well in-vitro, not necessarily those most relevant in-vivo.
- It provides no information about diversity or ecological relationships.

Targeted PCR / qPCR

Highly sensitive – but only for selected targets.

PCR-based testing is extremely useful when there is a specific clinical hypothesis (e.g., presence of potential pathogens such as *C. difficile*, *Giardia*).²

However:²

- You can only detect what you test for.
- It does not provide ecosystem-level composition.
- It does not assess diversity or resilience.
- It does not evaluate functional capacity.
- It cannot identify unexpected or emerging organisms.

PCR answers specific questions well, however it is not designed to assess ecosystem balance.

16S rRNA Gene Sequencing

Broader view, but limited resolution and no direct functional insight. 16S rRNA gene sequencing (16S) allows a broader survey of bacterial communities without culture. However, 16S rRNA sequencing is an amplicon-based method. It targets and amplifies a small region of the 16S ribosomal RNA gene, which is only present in bacteria and Archaea.

This creates several clinically relevant limitations:

- Often limited to genus-level resolution, with species-level identification carrying a high rate of false positives^{3 4}
- Inability to reliably distinguish clinically important species within the same genus.
- No detection of fungi or protists and 16S rRNA targets often miss Archaea^{3 4}
- It cannot identify novel or uncharacterised organisms
- No direct measurement of microbial functional genes or pathways.

Research directly comparing the two approaches confirms that 16S detects only part of the gut microbiota community revealed by shotgun sequencing, and that less abundant – but biologically meaningful, taxa are largely invisible to 16S methods.⁴ 16S can reveal broad trends, but it does not provide the resolution or functional depth required for high-confidence clinical interpretation.

Shotgun Metagenomics

Whole-community, species-level, function-aware measurement.
Shotgun metagenomics sequences all DNA in a sample.
This enables clinicians to assess:

- Composition: species and even strain-level identification of microorganisms
- Breadth: bacteria, fungi, Archaea, protists and other microbial DNA
- Discovery: detection of previously uncharacterised organisms
- Function: microbial genes and metabolic pathways

Studies demonstrate that shotgun metagenomic whole genome sequencing has multiple advantages over 16S amplicon methods, including enhanced detection of bacterial species, increased detection of diversity, and direct assessment of increased prediction functional of genes.^{5 6}



Comparison of microbiome testing methods

The table below summarises how common microbiome testing approaches differ in scope, resolution, and suitability for assessing ecosystem-level balance.

	CULTURE	TARGETED PCR/QPCR	16S RRNA GENES SEQUENCING	SHOTGUN METAGENOMICS
Detects Specific Pathogens	Yes (limited)	Yes (predefined targets)	No	Possible
Ecosystem-Level Composition	No	No	Partial (genus-level)	Yes (species-level)
Functional Insight	Partial (antimicrobial resistance only)	No	No (No direct functional data)	Yes (genes + pathways)
Suitable for Assessing Balance?	No	No	Limited	Yes

Ecosystem-level assessment reveals functional capacity

In an ecological model, clinical impact is determined not just by which organisms are present, but by what the ecosystem can do. This is called the Functional capacity.

These functions are distributed across multiple organisms that share metabolic pathways.



Understanding Functional Capacity

Functional capacity refers to the collective ability of the microbial community to carry out key biological activities.

Clinically relevant examples include

Short-chain fatty acid (SCFA) production (e.g. butyrate)

SCFAs are produced through bacterial fermentation of dietary fibre and play key roles in gut barrier integrity, immune modulation, and inflammation control. Reduced SCFA-producing capacity has been linked to increased intestinal permeability and systemic inflammatory burden.^{7 8}

Pro-inflammatory LPS signalling potential

LPS produced by Gram-negative bacteria triggers innate immune activation. Dysbiosis, with reduced diversity and overgrowth of pro-inflammatory LPS-producing organisms, shifts the microbiome towards a pro-inflammatory state, contributing to low-grade systemic inflammation and chronic disease risk.^{9 10}

Mucin degradation activity

Mucin-degrading organisms including *Akkermansia muciniphila* and *Bacteroides* species can promote mucus production when at optimal levels but can deplete the mucus layer at high levels. Excessive mucus degradation can compromise mucosal integrity and increase immune activation.^{11 12}

Protein fermentation by-products

Fermentation of undigested protein in the distal colon produces metabolites including trimethylamine, hydrogen sulphide, branched-chain amino acids, ammonia and phenolic compounds, some of which drive intestinal inflammation and may contribute to systemic toxic load at elevated concentrations.



Without functional information

Clinicians are left interpreting composition alone which often fails to explain chronic symptoms, relapse cycles, or systemic inflammatory burden.

Whole–ecosystem testing supports more informed clinical decision making

When clinicians move from targeted detection to whole–community measurement, the clinical questions evolve.


Instead of asking

- “Is a pathogen present?”
- “Is a specific organism elevated?”

Clinicians can now ask

- “Is this organism ecologically dominant or controlled?”
- “Is the ecosystem resilient or unstable?”
- “Is beneficial functional capacity reduced?”
- “Is inflammatory signalling potential elevated?”
- “How might diet, inflammation, or medication exposure be shaping this ecology?”

These are the questions that matter in complex, chronic presentations. Research linking gut dysbiosis to over 117 gastrointestinal and extra–gastrointestinal diseases underscore why understanding the ecosystem state, not just species presence, is clinically essential.¹³



Clinical-grade testing requires more than technology

Whole-community sequencing is necessary, but not sufficient.


For clinical application to be reliable and useful, it must be delivered through validated, accredited systems.

Metagenomics is only clinically meaningful when results are accurate, reproducible, and interpretable.

Technology alone doesn't deliver clarity. Integrated systems do. An integrated systems requires:

Validated sample collection and preservation

The microbiome changes rapidly once a sample is produced. Without robust preservation, results can be distorted by microbial growth, degradation, or DNA fragmentation.



Accredited laboratory standards (e.g., ISO 15189)

Clinical-grade testing requires operation to medical laboratory standards, including quality systems, traceability, and analytical validation.




High-quality bioinformatics pipelines

Metagenomic data is only as useful as the computational system interpreting it. Low-quality bioinformatics increases false positives, misclassification, and missing data.



Evidence-based interpretation frameworks

Results must be aligned with reproducible evidence in humans and presented in a clinically usable format.



Metagenomics is only clinically meaningful when results are accurate, reproducible, and interpretable.

A Practical Summary for Clinicians

Microbiome balance is an ecosystem property. It cannot be accurately assessed through partial or targeted testing alone.

Culture, PCR, and 16S methods can answer specific questions but they do not provide whole-ecosystem visibility. If the clinical goal is ecosystem-level understanding, whole-community metagenomics is required. And if the goal is actionable interpretation, the metagenomic assessments must be delivered through validated, clinical-grade systems.



REFERENCES

1. Rinninella E, et al. Gut microbiota, intestinal permeability, and systemic inflammation: a narrative review. *Internal and Emergency Medicine*. 2023. <https://link.springer.com/article/10.1007/s11739-023-03374-w>
2. Hilton SK, et al. Metataxonomic and metagenomic approaches vs. culture-based techniques for clinical pathology. *Frontiers in Microbiology*. 2016;7:484.
3. Rubiola S, et al. Comparison between 16S rRNA and shotgun sequencing in colorectal cancer, advanced colorectal lesions, and healthy human gut microbiota. *BMC Genomics*. 2024. <https://link.springer.com/article/10.1186/s12864-024-10621-7>
4. Comparri S, et al. Comparison between 16S rRNA and shotgun sequencing data for the taxonomic characterisation of the gut microbiota. *PubMed*. 2021. <https://pubmed.ncbi.nlm.nih.gov/33542369/>
5. Aggarwala V, et al. Analysis of the microbiome: Advantages of whole genome shotgun versus 16S amplicon sequencing. *PMC*. <https://pmc.ncbi.nlm.nih.gov/articles/PMC4830092/>
6. Aggarwala V, et al. *ibid.*
7. Portincasa P, et al. Short-Chain Fatty-Acid-Producing Bacteria: Key Components of the Human Gut Microbiota. *PMC*. 2023. <https://pmc.ncbi.nlm.nih.gov/articles/PMC10180739/>
8. Deleu S, et al. Short chain fatty acids and its producing organisms: An overlooked therapy for IBD? *eBioMedicine*. 2021. [https://www.thelancet.com/journals/ebiom/article/PIIS2352-3964\(21\)00086-4/fulltext](https://www.thelancet.com/journals/ebiom/article/PIIS2352-3964(21)00086-4/fulltext)
9. Candelli M, et al. Interaction between Lipopolysaccharide and Gut Microbiota in Inflammatory Bowel Diseases. *PMC*. 2021. <https://pmc.ncbi.nlm.nih.gov/articles/PMC8226948/>
10. Mohr AE, et al. Lipopolysaccharide and the gut microbiota: considering structural variation. *FEBS Letters*. 2022. <https://febs.onlinelibrary.wiley.com/doi/10.1002/1873-3468.14328>
11. Altamimi M, et al. Akkermansia muciniphila: A key player in gut microbiota-based disease modulation. *ScienceDirect*. 2025. <https://www.sciencedirect.com/science/article/pii/S0944501325002769>
12. Crost EH, et al. Breaking down barriers: is intestinal mucus degradation by Akkermansia muciniphila beneficial or harmful? *Infection and Immunity*. 2024. <https://journals.asm.org/doi/10.1128/iai.00503-24>
13. Taghinezhad-S S, et al. Understanding dysbiosis and resilience in the human gut microbio

Better Science.

Better Insights.

Better Health.

microba.com

AUSTRALIA
[1300 974 621](tel:1300974621)
hello@microba.com

UNITED KINGDOM
[+44 808 501 2748](tel:+448085012748)
hello.uk@microba.com

