

Overview: Measuring the microbiome as an ecosystem

The microbiome isn't a list of species. It's a living ecosystem. It should not be measured by testing fragments in isolation. The whole picture matters



A healthy gut microbiome requires a balanced ecosystem

Microbiome balance is not the property of a single microbe. It reflects the overall state of microbes in the gut ecosystem, shaped by:



Most testing technologies measure only fragments of this ecosystem, evaluating a partial list of well-known species and attempt to infer "balance" from that partial view. If the clinical question is ecosystem balance, the measurement must reflect the ecosystem.

Species-level detection

Measures organism presence in isolation — no diversity data, functional insight, or community behaviour

Ecosystem assessment

Captures diversity, abundance, functional capacity, and the physiological conditions shaping microbiome function.

Ecosystem-level assessment unlocks critical patient insights

In practice, microbiome testing is rarely about detection alone. The goal is interpretation: understanding how microbial ecology may be influencing:



Partial testing can't reliably answer critical clinical questions

- How does this organism's abundance compare to the overall ecosystem?
- Is the ecosystem resilient or vulnerable to instability?²³
- Is functional capacity reduced?
- Is it contributing to an anti-inflammatory or pro-inflammatory state?

Testing technologies vary in what they can measure and reveal

Culture-Based Testing



Useful for pathogen detection - not ecosystem assessment.

Culture can identify certain infectious organisms. Identifies organisms that grow in lab conditions
Misses many microbes that cannot be cultured
Provides no insight into diversity or function¹

Targeted PCR / qPCR



Highly sensitive, but only for selected targets.

- Excellent at confirming specific organisms
- Cannot assess broader community balance
- Does not detect organisms outside the selected panel.²

16S rRNA gene sequencing



Broader view but limited resolution and no direct functional insight.


- Surveys bacterial communities without culture
- Often limited to genus-level resolution
- Does not directly measure functional genes or pathways⁴

Shotgun Metagenomics



Whole-community, species-level, function-aware measurement.

- Sequences all DNA in the sample
- Detects bacteria, fungi, archaea and other microbes
- Provides species-level identification
- Enables functional gene and pathway profiling⁵



THE CRITICAL SHIFT:

Metagenomics captures the full genetic content of the microbial community, enabling ecosystem-level assessment.⁵ Not just selected fragments

Ecosystem-level assessment reveals functional capacity

Clinical impact is determined not just by which organisms are present, but by what the ecosystem produces. This is functional capacity. Without functional capacity insight, clinicians are left interpreting microbial composition alone, which often fails to explain chronic symptoms or relapsing presentation.

What is functional capacity ?

It is the collective ability of the gut microbiome to carry out key biological activities.⁵

Whole-ecosystem testing supports more informed clinical decision-making


Whole-ecosystem testing changes the questions clinicians can ask. Not just whether an organism is present – but whether the ecosystem is balanced, whether functional capacity is reduced, and whether it is showing signs of inflammation

Clinical-grade testing requires more than technology

Whole-community sequencing is necessary, but not sufficient.

For clinical application to be reliable and useful, it must be delivered through validated, accredited systems.

Clinical-grade testing requires:



Metagenomics is only clinically meaningful when results are accurate, reproducible, and interpretable. Technology alone doesn't deliver clarity. Integrated systems do.


The most accurate picture requires the right combination of technologies

No single method answers every clinical question. To deliver true ecosystem insight, testing must integrate:

Shotgun metagenomics for full community and functional capacity⁵

PCR for targeted pathogen confirmation

Validated Gut health markers to assess gut barrier function and inflammation



GI HEALTH MARKERS

TARGETED PCR + 16S SEQUENCING

SHOTGUN METAGENOMICS

Practical summary

Microbiome balance is an ecosystem property within the gut environment.¹²

Partial testing measures fragments, not systems. Functional capacity determines clinical impact.⁵ Ecosystem-level understanding requires whole-community measurement.

Clinical-grade clarity requires validated systems, not just sequencing. If the goal is confident, actionable insights, then the whole ecosystem must be.



Better Science.

Better Insights.

Better Health.

microba.com

AUSTRALIA
[1300 974 621](tel:1300974621)
hello@microba.com

UNITED KINGDOM
[+44 808 501 2748](tel:+448085012748)
hello.uk@microba.com

